




PRODUCT CATALOG

STAINLESS STEEL LABORATORY CASEWORK






Founded over 40 years ago, Inter Dyne Systems Inc. has become an innovator of stainless steel and acrylic laboratory furnishings. Fabricating stainless steel drying racks for laboratory glassware was only the beginning. Over time we gained a reputation as one of the leaders in the stainless steel laboratory accessories marketplace. We have built upon that foundation by continuing to offer high quality laboratory products with dependable service, consistent quality and on time delivery. Some of our offerings include but are not limited to stainless steel shelving and shelving components, stainless steel casework, countertops and tables, as well as innovative acrylic laboratory products such as dividers and splash guards.



4	Cabinet Specifications
4	Configuration Guide
5	Options Menu
6	Base Cabinet - Open Face
6	Base Cabinet - Door
7	Base Cabinet - Drawer
7	Base Cabinet - Door and Drawer
8	Base Cabinet - Corner
9	Sink Base Cabinet
9	Wall Mount Sink Enclosure
10	Mobile Cabinets
11	Hanging Cabinets
12	Tall Floor Cabinets
13	Wall Cabinets
14	Wall Cabinet - Corner
15	Filler Panels
16	Finishing Panel - Knee Space Apron
16	Finishing Panel - Knee Space Panel
16	Finishing Panel - Cabinet Backs
17	Accessories - Corner Sloped Top
17	Accessories - Sloped Top with End Caps
17	Accessories - Countertop Wall Brackets
18	Utility Table Leg Base
19	Utility Table Top
20	Scullery Sinks
21	Canopy Hood
21	Umbilicals
22	Stainless Steel Countertops



GAUGES OF BASE CABINET COMPONENTS

Aprons	18 ga.
Back Panels	18 ga.
Bottom Panels	18 ga.
Drawer/Door Outer Panel	20 ga.
Door Inner Panel	20 ga.
Drawer Bodies	20 ga.
Legs - 2" Sq. Tube	16 ga.
Shelves	18 ga.
Side Panels	18 ga.
Table Frames	18 ga.
Security Panels	20 ga.

CONFIGURATION GUIDE

SBC	1 or 2 Solid Swing Doors
SBC####OF	Open Face
SBC####4DW	4 Drawer Unit
SBC####A	1 Door, 1 Drawer or 2 Door, 2 Drawer
SBC####B	1 Drawer, 2 Door
-ADA	ADA Height 31.750" Tall
-SK	Sink Base Cabinet
-CC	Corner Cabinet
-GD	Framed Glass Door
-RH	Right Hinge
-LH	Left Hinge

All surfaces shall be stainless steel. Specific test data is available upon request. Type 304 with a #4 finish standard.

For complete specifications, request our full line of product specification. Subject to change.

DOOR & DRAWER FRONTS

Standard Default Inset Steel
(Optional Overlay)



DRAWER SLIDE

Standard Default 100 lb. Full
Extension, All Stainless Steel



MAGNETIC CATCH

Standard Default Magnetic Catch
(Optional Roller Catch)



PULL

Stainless Steel
4" Wire Pull
(Optional Recessed Pull)




HINGE


Inset Steel
Stainless Steel 5 Knuckle



GLASS FRAMED DOOR

- Standard 1/4" Clear Float Glass
- Optional 1/4" Laminated Safety Glass
- Optional 1/4" Tempered Safety Glass





Standing height base cabinets are 34.75" high and 22.375" deep. Units are available in various standard widths as shown. Door units feature a single shelf as standard. Additional shelves are available. Specify left- or right-hinged doors using the suffix "LH" or "RH". Finished back panels are available.

OPEN FACE
BASE CABINET



WIDTH	ITEM #
18"	SBC3518-OF
24"	SBC3524-OF
30"	SBC3530-OF
36"	SBC3536-OF
42"	SBC3542-OF
48"	SBC3548-OF

OPEN FACE
BASE CABINET - ADA

WIDTH	ITEM #
18"	SBC3218-OF-ADA
24"	SBC3224-OF-ADA
30"	SBC3230-OF-ADA
36"	SBC3236-OF-ADA
42"	SBC3242-OF-ADA
48"	SBC3248-OF-ADA

DOOR BASE CABINET




WIDTH	ITEM #
18"	SBC3518-LH
18"	SBC3518-RH
24"	SBC3524-LH
24"	SBC3524-RH
30"	SBC3530
36"	SBC3536
42"	SBC3542
48"	SBC3548

DOOR BASE CABINET - ADA

WIDTH	ITEM #
18"	SBC3218-LH-ADA
18"	SBC3218-RH-ADA
24"	SBC3224-LH-ADA
24"	SBC3224-RH-ADA
30"	SBC3230-ADA
36"	SBC3236-ADA
42"	SBC3242-ADA
48"	SBC3248-ADA

Standing height base cabinets are 34.75" high, ADA Height of 31.75" and 22.375" deep. Units are available in various widths as shown. Drawer units feature 100 lb. full extension stain-less steel slides as standard. Door units feature a single shelf as standard. Additional shelves are available. Specify left- or right-hinged doors using the suffix "LH" or "RH". If locks are required use the suffix "LK". Finished back panels are available.



DRAWER BASE
CABINET



WIDTH	ITEM #
18"	SBC3518-4DW
24"	SBC3524-4DW
30"	SBC3530-4DW
36"	SBC3536-4DW
42"	SBC3542-4DW
48"	SBC3548-4DW

DRAWER BASE
CABINET - ADA

WIDTH	ITEM #
18"	SBC3218-4DW-ADA
24"	SBC3224-4DW-ADA
30"	SBC3230-4DW-ADA
36"	SBC3236-4DW-ADA
42"	SBC3242-4DW-ADA
48"	SBC3248-4DW-ADA


DOOR AND DRAWER
BASE CABINET



WIDTH	ITEM #
18"	SBC3518-A-LH
18"	SBC3518-A-RH
24"	SBC3524-A-LH
24"	SBC3524-A-RH
30"	SBC3530-A
36"	SBC3536-A
42"	SBC3542-A
48"	SBC3548-A

DOOR AND DRAWER
BASE CABINET - ADA

WIDTH	ITEM #
18"	SBC3218-A-LH-ADA
18"	SBC3218-A-RH-ADA
24"	SBC3224-A-LH-ADA
24"	SBC3224-A-RH-ADA
30"	SBC3230-A-ADA
36"	SBC3236-A-ADA
42"	SBC3242-A-ADA
48"	SBC3248-A-ADA



Standing height base cabinets are 34.75" high and 22.375" deep. Units are available in various standard widths as shown. Drawer units feature 100 lb. full extension stainless steel slides as standard. Door units feature a single shelf as standard. Additional shelves are available. Specify left- or right-hinged doors using the suffix "LH" or "RH". Finished back panels are available.

BASE CABINET - 1 DRAWER
2 DOOR

WIDTH	ITEM #
30"	SBC3530-B
36"	SBC3536-B
42"	SBC3542-B
48"	SBC3548-B




BASE CABINET - 1 DRAWER
2 DOOR - ADA

WIDTH	ITEM #
30"	SBC3230-B-ADA
36"	SBC3236-B-ADA
42"	SBC3242-B-ADA
48"	SBC3248-B-ADA

CORNER BASE CABINET




ITEM #
SBC3514-CC-LH
SBC3514-CC-RH
SBC3530-4DW
SBC3536-4DW
SBC3542-4DW
SBC3548-4DW



SINK BASE CABINET


Standing height base cabinets are 34.75" high and 22.375" deep. Units are available in various widths as shown. Specify left- or right-hinged doors using the suffix "LH" or "RH". Back side is open to allow for plumbing.

WIDTH	ITEM #
24"	SBC3524-SK-LH
24"	SBC3524-SK-RH
30"	SBC3530-SK
36"	SBC3536-SK
42"	SBC3542-SK
48"	SBC3548-SK
60"	SBC3560-SK
48"	SBC3548-SKI
60"	SBC3560-SKI



SINK BASE CABINET - ADA

ADA sink base cabinets are 31.75" high and 22.375" deep. Units are available in various widths as shown.



WIDTH	ITEM #
36"	SBC3236-SKA
42"	SBC3242-SKA
48"	SBC3248-SKA
30"	SBP3236
36"	SBP3242
42"	SBP3248
48"	SBP3260

WALL MOUNT SINK ENCLOSURE - ADA

WIDTH-DEPTH	ITEM #
30"-23"	ADSE2230-23
30"-24"	ADSE2230-24
30"-29"	ADSE2230-29
36"-23"	ADSE2236-23
36"-24"	ADSE2236-24
36"-29"	ADSE2236-29
48"-23"	ADSE2248-23
48"-24"	ADSE2248-24
48"-29"	ADSE2248-29

Stainless steel sink enclosures are supplied with 12-gauge end support brackets, an 18-gauge apron, and 18-gauge removable front panel.





Mobile cabinets are available in a sitting 28.75" height or a 31.75" ADA height. They come standard with a finished back and finished top. They come standard with a 3" locking caster. All mobile cabinets with drawers come standard with a tip resistant counter weight. Low profile and standing height available on request.

MOBILE CABINET DOOR UNIT
ITEM #
SMC2918-LH
SMC2918-RH
SMC2924-LH
SMC2924-RH
SMC2930
SMC2936
SMC2942
SMC2948

MOBILE CABINET - ADA DOOR UNIT
ITEM #
SMC3218-LH-ADA
SMC3218-RH-ADA
SMC3224-LH-ADA
SMC3224-RH-ADA
SMC3230-ADA
SMC3236-ADA
SMC3242-ADA
SMC3248-ADA

MOBILE CABINET DOOR/DRAWER UNIT
SMC2918-A-LH
SMC2918-A-RH
SMC2924-A-LH
SMC2924-A-RH
SMC2930-A
SMC2936-A
SMC2942-A
SMC2948-A

MOBILE CABINET - ADA DOOR/DRAWER UNIT
SMC3218-A-LH-ADA
SMC3218-A-RH-ADA
SMC3224-A-LH-ADA
SMC3224-A-RH-ADA
SMC3230-A-ADA
SMC3236-A-ADA
SMC3242-A-ADA
SMC3248-A-ADA

MOBILE CABINET 3 DRAWER UNIT
SMC2918-3DW
SMC2924-3DW
SMC2930-3DW
SMC2936-3DW
SMC2942-3DW
SMC2948-3DW

MOBILE CABINET - ADA 3 DRAWER UNIT
SMC3218-3DW-ADA
SMC3224-3DW-ADA
SMC3230-3DW-ADA
SMC3236-3DW-ADA
SMC3242-3DW-ADA
SMC3248-3DW-ADA

Hanging cabinets are available in a height of 24.875" or ADA height of 27.875". Finished backs are optional.



HANGING CABINET DOOR UNIT
ITEM #
SHC2518-LH
SHC2518-RH
SHC2524-LH
SHC2524-RH
SHC2530
SHC2536
SHC2542
SHC2548

HANGING CABINET - ADA DOOR UNIT
ITEM #
SHC2818-LH-ADA
SHC2818-RH-ADA
SHC2824-LH-ADA
SHC2824-RH-ADA
SHC2830-ADA
SHC2836-ADA
SHC2842-ADA
SHC2848-ADA

HANGING CABINET DOOR/DRAWER UNIT
SHC2518-A-LH
SHC2518-A-RH
SHC2524-A-LH
SHC2524-A-RH
SHC2530-A
SHC2536-A
SHC2542-A
SHC2548-A

HANGING CABINET - ADA DOOR/DRAWER UNIT
SHC2818-A-LH-ADA
SHC2818-A-RH-ADA
SHC2824-A-LH-ADA
SHC2824-A-RH-ADA
SHC2830-A-ADA
SHC2836-A-ADA
SHC2842-A-ADA
SHC2848-A-ADA

HANGING CABINET 3 DRAWER UNIT
SHC2518-3DW
SHC2524-3DW
SHC2530-3DW
SHC2536-3DW
SHC2542-3DW
SHC2548-3DW

HANGING CABINET - ADA 3 DRAWER UNIT
SHC2818-3DW-ADA
SHC2824-3DW-ADA
SHC2830-3DW-ADA
SHC2836-3DW-ADA
SHC2842-3DW-ADA
SHC2848-3DW-ADA



Tall floor cabinets are 84" high and available in a standard depth of 22.375". Units are available in various widths as shown. Door units feature five shelves (four adjustable and one fixed) as standard. Additional shelves are available. Specify left- or right-hinged doors using the suffix "LH" or "RH" after the two digit depth number for cabinets under 30" width.

OPEN FACE FLOOR CABINET



WIDTH	ITEM #
18"	SFC8418-OF
24"	SFC8424-OF
30"	SFC8430-OF
36"	SFC8436-OF
42"	SFC8442-OF
48"	SFC8448-OF

DOOR FLOOR CABINET

WIDTH	ITEM #
18"	SFC8418-LH
18"	SFC8418-RH
24"	SFC8424-LH
24"	SFC8424-RH
30"	SFC8430
36"	SFC8436
42"	SFC8442
48"	SFC8448




FRAMED GLASS DOOR FLOOR CABINET




WIDTH	ITEM #
18"	SFC8418-GD-LH
18"	SFC8418-GD-RH
24"	SFC8424-GD-LH
24"	SFC8424-GD-RH
30"	SFC8430-GD
36"	SFC8436-GD
42"	SFC8442-GD
48"	SFC8448-GD

Wall cabinets are available in 24", 30" and 36" heights and come standard with a 13" depth. For 16" deep cabinets, use the suffix "16". Units are available in various widths as shown. They feature two shelves as standard for 30" units. Specify left- or right-hinged doors using the suffix "LH" or "RH" after the two digit depth number.



OPEN FACE WALL CABINET

WIDTH	24" ITEM #	30" ITEM #	36" ITEM #
18"	SWC2418-OF	SWC3018-OF	SWC3618-OF
24"	SWC2424-OF	SWC3024-OF	SWC3624-OF
30"	SWC2430-OF	SWC3030-OF	SWC3630-OF
36"	SWC2436-OF	SWC3036-OF	SWC3636-OF
42"	SWC2442-OF	SWC3042-OF	SWC3642-OF
48"	SWC2448-OF	SWC3048-OF	SWC3648-OF



DOOR WALL CABINET



WIDTH	24" ITEM #	30" ITEM #	36" ITEM #
18"	SWC2418-LH	SWC3018-LH	SWC3618-LH
18"	SWC2418-RH	SWC3018-RH	SWC3618-RH
24"	SWC2424-LH	SWC3024-LH	SWC3624-LH
24"	SWC2424-RH	SWC3024-RH	SWC3624-RH
30"	SWC2430	SWC3030	SWC3630
36"	SWC2436	SWC3036	SWC3636
42"	SWC2442	SWC3042	SWC3642
48"	SWC2448	SWC3048	SWC3648

FRAMED GLASS DOOR WALL CABINET

WIDTH	24" ITEM #	30" ITEM #	36" ITEM #
18"	SWC2418-GD-LH	SWC3018-GD-LH	SWC3618-GD-LH
18"	SWC2418-GD-RH	SWC3018-GD-RH	SWC3618-GD-RH
24"	SWC2424-GD-LH	SWC3024-GD-LH	SWC3624-GD-LH
24"	SWC2424-GD-RH	SWC3024-GD-RH	SWC3624-GD-RH
30"	SWC2430-GD	SWC3030-GD	SWC3630-GD
36"	SWC2436-GD	SWC3036-GD	SWC3636-GD
42"	SWC2442-GD	SWC3042-GD	SWC3642-GD
48"	SWC2448-GD	SWC3048-GD	SWC3648-GD




CORNER WALL CABINET OPEN FACE
ITEM #
SWC3018-OF-CC



CORNER WALL CABINET DOOR
ITEM #
SWC3018-OF-CC-LH
SWC3018-OF-CC-RH



CORNER WALL CABINET GLASS DOOR
ITEM #
SWC3018-GD-CC-LH
SWC3018-GD-CC-RH



Base cabinet filler panels are available in 32" and 35" height. Floor case filler panels are available in 84" height. They are all available in any width required by entering the dimension in lieu of the "XX" in the part number. Note: Inter Dyne base cabinets are 22.375" deep.



End Fill

ITEM #
SEFP-0132
SEFP-0135




Base Cabinet Front Fill

ITEM #
SBCFP-0132
SBCFP-0135



Rear Fill

ITEM #
SRFP-0132
SRFP-0135



Base Corner Fill

ITEM #
SBCFP-0132-C
SBCFP-0135-C



Floor Case Front Fill

ITEM #
SFCFFP-0132-C
SFCFFP-0135-C



Floor Case Rear Fill

ITEM #
SFCRFP-0132-C
SFCRFP-0135-C



Wall Cabinet Front Fill

ITEM #
SWCFP-24
SWCFP-30
SWCFP-36



Wall Cabinet Corner Fill

ITEM #
SWCFP-24-C
SWCFP-30-C
SWCFP-36-C

KNEE SPACE APRONS AND KNEE SPACE PANELS

ITEM #	ITEM #	ITEM #
KSA-18	SKP-30	SKP-30-B
KSA-24	SKP-36	SKP-36-B
KSA-30	SKP-42	SKP-42-B
KSA-36	SKP-48	SKP-48-B
KSA-42	SKP-54	SKP-54-B
KSA-48	SKP-60	SKP-60-B
KSA-54		
KSA-60		




KNEE SPACE APRON



KNEE PANEL

-B BOTH

CABINET FINISHED BACKS



Island finished back panels are available upon request. Call for part number and cost. Part number example is FB-SBC____ (width then depth) 18 gauge panel.

CORNER SLOPED TOP



ITEM #
SST-CORNER13-SS
SST-CORNER16-SS

SLOPED TOP WITH END CAPS


ITEM #
SST13-18
SST13-24
SST13-30
SST13-36
SST13-42
SST13-48
SST16-18
SST16-24
SST16-30
SST16-36
SST16-42
SST16-48

Sloped tops are 6" tall unless otherwise noted. Sloped tops can attach to cabinet tops or rear wall.



COUNTERTOP WALL BRACKETS

12 gauge stainless steel wall brackets for countertops. Brackets are 18" high and 22" or 28" deep. Supplied with installation holes. Part number reflects non-installation side (example: CTB-24-LH, would have install flanges on the right side).



ITEM #	DEPTH	HEIGHT
CTB-24-LH	22"	18"
CTB-24-RH	22"	18"
CTB-30-LH	28"	18"
CTB-30-RH	28"	18"


All stainless steel utility tables are made with 2" square tube legs and adjustable bullet feet. Front to back and center stretchers with fully welded 4" tall aprons are standard. Center stretchers can be added to the rear at no extra cost. Adjustable legs are available and fixed heights are 35" tall, unless otherwise noted. Drawer options and lockable casters are available upon request. Adjustable height tables are manually adjustable with supplied locking clevis pins. *Add "CAST" to Item # for four locking casters.

UTILITY TABLE LEG BASE				
	FIXED ITEM #	ADJUSTABLE ITEM #	LENGTH	DEPTH
	STF-3024-36	STF-3024-ADJ	28"	22"
	STF-3030-36	STF-3030-ADJ	28"	28"
	STF-3036-36	STF-3036-ADJ	28"	34"
	STF-3624-36	STF-3624-ADJ	34"	22"
	STF-3630-36	STF-3630-ADJ	34"	28"
	STF-3636-36	STF-3636-ADJ	34"	34"
	STF-4824-36	STF-4824-ADJ	46"	22"
	STF-4830-36	STF-4830-ADJ	46"	28"
	STF-4836-36	STF-4836-ADJ	46"	34"
	STF-5424-36	STF-5424-ADJ	52"	22"
	STF-5430-36	STF-5430-ADJ	52"	28"
	STF-5436-36	STF-5436-ADJ	52"	34"
	STF-6024-36	STF-6024-ADJ	58"	22"
	STF-6030-36	STF-6030-ADJ	58"	28"
	STF-6036-36	STF-6036-ADJ	58"	34"
	STF-6624-36	STF-6624-ADJ	64"	22"
	STF-6630-36	STF-6630-ADJ	64"	28"
	STF-6636-36	STF-6636-ADJ	64"	34"
	STF-7224-36	STF-7224-ADJ	70"	22"
	STF-7230-36	STF-7230-ADJ	70"	28"
	STF-7236-36	STF-7236-ADJ	70"	34"
	STF-8424-36	STF-8424-ADJ	82"	22"
	STF-8430-36	STF-8430-ADJ	82"	28"
	STF-8436-36	STF-8436-ADJ	82"	34"
	STF-9624-36	STF-9624-ADJ	94"	22"
	STF-9630-36	STF-9630-ADJ	94"	28"
	STF-9636-36	STF-9636-ADJ	94"	34"

Stainless steel table tops are 1" thick 16-gauge with stainless steel hatchanel substrate and sprayed with a sound deadening compound. Table tops are available in 24", 30", and 36" depths as standard. Customizations are available.

UTILITY TABLE TOPS


FIXED ITEM #	LENGTH	DEPTH	WIDTH
STT-3024-1	30"	24"	1"
STT-3624-1	36"	24"	1"
STT-4224-1	42"	24"	1"
STT-4824-1	48"	24"	1"
STT-5424-1	54"	24"	1"
STT-6024-1	60"	24"	1"
STT-6624-1	66"	24"	1"
STT-7224-1	72"	24"	1"
STT-8424-1	84"	24"	1"
STT-9624-1	96"	24"	1"
STT-3030-1	30"	30"	1"
STT-3630-1	36"	30"	1"
STT-4230-1	42"	30"	1"
STT-4830-1	48"	30"	1"
STT-5430-1	54"	30"	1"
STT-6030-1	60"	30"	1"
STT-6630-1	66"	30"	1"
STT-7230-1	72"	30"	1"
STT-8430-1	84"	30"	1"
STT-9630-1	96"	30"	1"
STT-3036-1	30"	36"	1"
STT-3636-1	36"	36"	1"
STT-4236-1	42"	36"	1"
STT-4836-1	48"	36"	1"
STT-5436-1	54"	36"	1"
STT-6036-1	60"	36"	1"
STT-6636-1	66"	36"	1"
STT-7236-1	72"	36"	1"
STT-8436-1	84"	36"	1"
STT-9636-1	96"	36"	1"



Length

Depth

UNDERSIDE VIEW



EXPLODED VIEW

Stainless steel 14-gauge scullery sinks come with 2" square tube legs on ADJ bullet feet with front to back and center stretchers. Tops are pitched for drainage. The "RFD" style comes with a lower shelf and boxed curb. Umbilical and column notches are available. Any size or configuration is available. Comes with wall clips, drain outlets, and basket strainers.




ITEM #	LENGTH	DEPTH	# OF SINKS
SCUL-72-1B	72"	30"	1
SCUL-72-2B	72"	30"	2
SCUL-96-2B	96"	30"	2
SCUL-96-3B	96"	30"	3

ITEM #	LENGTH	DEPTH	# OF SINKS
RFDS-72-18	72"	33"	1
RFDS-72-28	72"	33"	2
RFDS-96-28	96"	33"	2
RFDS-96-38	96"	33"	3
RFDS-120-28	120"	33"	2
RFDS-120-38	120"	33"	3



STAINLESS STEEL SCULLERY SPECS




- 1B - 1 Bowl
- 2B - 2 Bowls
- 3B - 3 Bowls

CONSTRUCTION


- 14 Gauge
- Type 304 Stainless Steel
- Polished to a #4 Finish

ADDITIONAL FEATURES


- Lower Shelf
- Center Set of Legs Near Sink
- Lever Handle Drain
- Tall Front Apron



PROVIDED WALL CLIP




MARINE EDGE



All stainless steel 16-gauge canopy hoods (wall mount or island style) are supplied with integral brackets to attach ceiling rods with duct collar and gutter. Standard and custom sizes available. Baffles available upon request.


ISLAND CANOPY HOOD

ITEM #	LENGTH	DEPTH	HEIGHT
ICH-3630-18	36"	30"	18"
ICH-4830-18	48"	30"	18"
ICH-6030-18	60"	30"	18"
ICH-7230-18	72"	30"	18"
ICH-9630-18	96"	30"	18"



WALL CANOPY HOOD

ITEM #	LENGTH	DEPTH	HEIGHT
WCH-3630-18	36"	30"	18"
WCH-4830-18	48"	30"	18"
WCH-6030-18	60"	30"	18"
WCH-7230-18	72"	30"	18"
WCH-9630-18	96"	30"	18"



Stainless steel utility chases "Umbilicals" are 18-gauge panels with unit strut interiors. Supplied with removable panels and mounting collars.

ISLAND UMBILICAL

ITEM #
SIU-2672
SIU-2684
SIU-2696



WALL UMBILICAL

ITEM #
SWU-2672
SWU-2684
SWU-2696




Stainless steel countertops and tables are constructed from type 304 stainless steel, polished to a #4 finish. Type 316 stainless steel is available upon request for an additional charge. Countertops and tables are constructed from 16-gauge thick material and sink bowls from 14-gauge material.

Table frames are made using 2" square tubular construction with adjustable bullet-type feet.

- 14-Gauge Sink Bowls
 - Drain Outlets With Basket Strainers
 - Optional Wall Mount Brackets
 - Wall Clips Provided
 - Pitches Available
 - Partial & Full End Splashers
- Tops Crated for Shipping
 - Price Includes A Submittal Drawing
 - Stainless Steel Hat Channel Substrate
 - Sound Deadening Compound Applied



STAINLESS STEEL CUSTOM COUNTERTOPS




MARINE EDGE

MECHANICAL JOINT

The Victoria or "V" pegboard is constructed of lightweight stainless steel, offered in a range of sizes and variety of mounting options. Along the bottom is a drip trough that funnels water away to a drain, allowing equipment to dry quickly and remain sanitized. Each "V" style pegboard includes a drip trough, wall bracket, standard 6" white pegs and 3 ft. of clear plastic drain hose.

The "B" Baron is also made of lightweight stainless steel. It is offered in three sizes and a variety of mounting options. A drip deflector at the bottom diverts water away from walls and directs it to a sink or drain. Each "B" style pegboard includes a wall bracket and standard 6" white pegs.




"V" VICTORIA
STYLE

"B" BARON
STYLE

"V" VICTORIA STYLE PEGBOARD



ITEM #	W x H	# PEGS	SHIP WT.
V1824	18" x 24"	15	13 lbs
V2418	24" x 18"	16	13 lbs
V2424	24" x 24"	20	16 lbs
V2430	24" x 30"	32	23 lbs
V2436	24" x 36"	40	25 lbs
V3024	30" x 24"	25	23 lbs
V3030	30" x 30"	50	24 lbs
V3036	30" x 36"	40	23 lbs
V3624	36" x 24"	30	25 lbs
V3630	36" x 30"	60	23 lbs
V3636	36" x 36"	66	25 lbs
V4824	48" x 24"	40	24 lbs
V4830	48" x 30"	48	28 lbs
V4836	48" x 36"	88	31 lbs



WALL MOUNT

FINISHED BACK

FREE STANDING ASSEMBLY

RAISED TOP ASSEMBLY

WALL STANDARD MOUNT



Prime industries, inc.

4611 Main St., Lisle, IL 60532

Office: (630) 725-9200 • Fax: (630) 725-9201 • sales@piilab.com

www.piilab.com