


Stainless Steel Carboy Shelves

Wall mount carboy shelf



<i>Features</i>	<i>Model Number</i>	<i>Dimensions</i>
• All stainless steel construction	CS-18-W	18"w x 12"d
• Made from 18 gauge material	CS-20-W	20"w x 12"d
• Type 304 stainless steel, #4 finish	CS-24-W	24"w x 12"d
	CS-30-W	30"w x 12"d
• 4"h wraparound	CS-32-W	32"w x 12"d
	CS-36-W	36"w x 12"d
	CS-48-W	48"w x 12"d



INTER DYNE
SYSTEMS

Inter Dyne Systems, Inc.
676 Ellis Road, Norton Shores, MI 49441
Ph: 877.4ss.labs
www.interdynesystems.com